

Hardwired Pad MONITOR

A-BME-03FAQ'S

Q. How do I connect the pad to the monitor?

R. Simply plug the cable that looks like a landline telephone plug into the monitor

Q. Can I use batteries and a power supply at the same time?

R. Yes, in fact it is always recommended to do so in case of a power outage. The batteries will act as a power backup

Q. How do I store the monitor when it is not in use?

R. Disconnecting the pad and the power supply is ideal. Also taking the batteries out is also important as when they start to go flat the alarm will continue to beep.

Q. Can this be used for aged care if the patient goes to respite or hospital

R. Yes. There are some interconnection cables required but this could operate as a part of a nurse call system.

Technical Specifications:

- Pause mode to attend to a patient without triggering alarm.
- Tamper proof reset option using magnetic caregiver key.
- Choice of tone, and volume control for alarm.
- 0, 2s, and 15min. alarm delay settings.
- Port to connect with nurse call system (NC/NO).
- Operates on batteries or via AC adaptor.
- Low battery notification.
- Input lost notification (pad lost,