

CURA1 TRANSMITTER A-CL433-TX FAQ'S

Q. Can I use the one transmitter for two different pads?

- R. Yes, The transmitter fits to any Cura1 Wireless pad or mat.
- Cordi How far will the transmitter signal? Q.

Y 11.11

The distance is between 20-50 meters in direct line of sight however we R. recommend to test the distance as other radio signals can sometimes ce compli interfere.

Rules. Oper Can you replace the battery when it gets low battery? Q.

Although the battery inside the compartment looks like a standard coin R. battery, the firmware has a non- tamper program so that the change of battery will not extend the life of the transmitter. device accept

lived, including ence, and (2) any interference Cause un Technical Specifications: nce that may

oper 2mm x1.2mm 4pin pitch plug

- 433MHz operating frequency
- Up to a 20m range*
- Uses a Lith-ION long life battery 1 year battery life 6 month warranty
- Will go into sleep mode when not in use to save battery life
- Status/in-use LED indication
- Communicates with A- CBM- bedside controller
- Incorporates on/off
- Fits within floor, bed & chair sensor transmitter enclosure

Smart Caller Safe Life 5/270 Lower Dandenong Rd Mordialloc Vic 3195 PH: +61 3 9588 0833 E: sales@safelife.com.au