

cura¹
Fall Prevention Products



A fully featured wander detection system ideal to notify when a resident known for Resident to Resident Abuse (RRA)* is within a pre-determined range of another specific resident's room.

The DoorAlert LT RRA Monitor is a plug and play system designed to assist carers in preventing wandering of (typically) dementia sufferers.

Automatically detects wandering person and alerts carers.

DoorAlert LT RRA Monitor

Features:

- *Automatically detects wandering person and alerts carers.
- *Fully wireless – no cabling required for easy installation, configuration, operation. Requires a power point for the Door Alert LT RRA monitor.
- *Fast, simple 'plug & play' setup and installation- no specialist installers required.
- *Triggers call point via your chosen cura¹ cordless monitor.
- *Door Alert LT RRA Monitor emits an audible alarm.
- *Robust, stable, and consistent performance.
- *Only the monitored person is detected. Alarms can not be caused by pets or other residents.
- *Alarm can notify when person is in alarmed area or only when the door is opened.
- *Detection range adjustable from 0-5m.
- *Flexibility: the system can be configured in a variety of ways and expanded to use additional alarm devices through cordless monitor.
- *Tamper-proof installation.
- *Transmitter low battery indication.
- *Transmitter can be worn on the wrist, as a pendant or attached to clothing.
- *Transmitter has guaranteed 12 months continuous operation.
- *Escort mode allows a carer to escort a monitored person through a door without activation and alarm if required.



SafeLife
sales@safelife.com.au
03-9588-0392

**Resident-to-Resident abuse or aggression (RRA) has been defined as 'negative and aggressive physical, sexual or verbal interactions between long-term care residents that in a community setting would likely be construed as unwelcome and high potential to cause physical or psychological distress in the recipient'.*


[Prevalence and nature of resident-to-resident abuse incidents in Australian residential aged care](#) [Catherine M. Joyce](#)

cura: GUARANTEE: All Products supplied with full warranty and guaranteed to be free from defects in material and workmanship. The range of wireless devices may vary and be dependent on your specific environmental conditions. cura: products are designed to augment well designed care programs and practices, cura: does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision and are at risk.




HOW TO ORDER:

The DoorAlert LT RRA, with your choice of wireless monitor to connect to nurse call systems, is all that is required to have a fully functional Resident to Resident Abuse (RRA)* wander system.

DoorAlert LT Kit:








Picture	Code	Description	Whats in the box
	4000RRA	DoorAlert LT RRA. Automatically detects wandering person and alerts caregivers. Flashing LED and adjustable audible alarm.	1 x Door Monitor, 1 x Universal Transmitter, 1 x Wristband Transmitter, 1 x Plastic Wrist Strap, 1 x White Buckle Wrist Strap, 1 x Lanyard W/Safety Break Clip, 1 x 12Vdc Plug Pack, 1 x Magnetic Door Reed Switch, 1 x Caregiver Reset Key, 1 x Quick Start Guide.


Cordless Monitors:

Picture	Code	Description
	2544 or 2549	Multiport Expander. Plugs directly into the socket of your call point and alerts your nurse call system when the DoorAlert LT RRA alarms. The rooms resident can be monitored with another 4 wireless devices if required including a wearable pendant. Silent within the room. 2549 version for different nurse call systems or hard to reach call points.
	2700	6 Channel Monitor. Connects to your nurse call system via a connection lead. Can have another 5 wireless monitoring devices programmed. Can be silent within the room or emits an audible alarm.
	2702	3 Channel Monitor. Connects to your nurse call system via a connection lead. Can have another 2 wireless monitoring devices programmed. Can be silent within the room or emits an audible alarm.

Please note: Choose other monitoring devices from the cordless Fall Prevention Range. Contact us for a copy of our cura¹ product line.

Replacement Accessories:

Picture	Code	Description
	4001	DoorAlert LT RRA Wristband Transmitter. Guaranteed 12 months continuous operation. Includes strap. Lightweight, small and showerproof to IP66. LED indication of low battery conditions.
	2569	Plastic Wrist Strap. Replacement wrist strap for residents / patient's wristband transmitters.
	2589	White Buckle Wrist Strap. Replacement wrist straps for DoorAlert LT RRA wristband transmitter.
	2642	Caregiver Reset Key. For use by caregivers to reset the Door Monitor alarms at the monitor. Economical and easy to use.
	2567	Lanyard With Safety Break Clip. Lanyard for DoorAlert LT RRA wristband transmitter.
	4002	DoorAlert LT RRA Reset Button. For use by caregivers to reset the Door Monitor alarms. Fully wireless wall mount installation. Up to 10m range. Shower proof IP66 rating. Not included in DoorAlert LT RRA.
	2404	Plug Pack 12Vdc. Suitable for DoorAlert LT RRA Monitor.



cura¹ GUARANTEE: All Products supplied with full warranty and guaranteed to be free from defects in material and workmanship. The range of wireless devices may vary and be dependent on your specific environmental conditions. cura¹ products are designed to augment well designed care programs and practices, cura¹ does not claim that any of its products are suitable as a substitute, or as a total solution for individuals that require care and constant supervision and are at risk.