


2669 Cordless Under Mattress **PreAlert** Bed Pad

cura¹
Helping To Care



*... alerts earlier so you can
respond faster*

2669 Cordless Under Mattress **PreAlert** Bed Pad

A caregiver has more time to prevent a potential fall when they receive an early warning of a patient or resident sitting or rolling to the edge of the bed.

- Can provide an early warning when a patient or resident sits on the edge of the bed.
- Installed under the mattress for maximum discretion & user comfort.
- Self-calibration technology, zero up the weight of the mattress.
- Works with a wide variety of mattresses.
- Full width single bed pad (810 x 240 mm) means reliable operation with practically no false alarms.
- The pad incorporates a delay function in the centre area to prevent false alarms when repositioning.
- Sitting on the edge to get into bed does NOT cause an alarm.
- Transmitter included.

*The best solution for an
early warning of out of bed alarm*



Default delay setting in the center
reduces the incidence of false alarms.



Repositioning in bed such as rolling
to one side will not cause an alarm.



Sensor pad initiates the alarm when
patient sits up and moves toward the
edge of the bed.



Please Note: Pad is installed
under the mattress but is
shown on top for
illustration purposes.

WHAT'S IN THE BOX (2669)



2669 Cordless Under
Mattress PreAlert Bed Pad

Convenient foldable
design reduces transport
& storage costs.



COMPATIBLE PRODUCTS

- 3508 LED Pager
- 2544/9 Cordless Multiport Expander
- And all other cura¹ cordless receivers.



Certifications may vary between products.

cura¹
Helping To Care

10 Gray Street, Kilkenny, South Australia 5009

sales@cura1.com | www.cura1.com

+61 (0) 8 8268 8343 (International)

1300 12-CARE (2273)

PB-025