

Hardwired Bed Pad

A-BPLB-01 FAQ'S

Q. How do I connect the pad to the monitor?

R. Simply plug the cable that looks like a landline telephone plug into the monitor

Q. Can I clean the Pad?

R. You cannot wash the pad where it is submerged in water, however you can wipe the pad with a cleaning solution and damp cloth.

Q. How do I store the pad when it is not in use?

R. Although you cannot fold or roll the pad, you can lay it flat. Protecting the sheet with a thick layer of cardboard is ideal.

Q. Can this be used for aged care if the patient goes to respite or hospital

R. Yes. There are some interconnection cables required but this could operate as a part of a nurse call system.

Technical Specifications:

- Simple to operate
- 760 x 510 mm
- 1 Year life expectancy
- Easy to clean maintenance
- For use with Cura1 BME
- weight range for reliable operation: 30 kg to 200 kg
- Made from soft PVC
- Vent holes which reduce the risk of moisture
- Replaceable Break cord

Smart Caller Safe Life
5/ 270 Lower Dandenong Rd Mordialloc Vic
3195
PH: +61 3 9588 0833
E: sales@safelife.com.au