

## Hardwired Bed, Chair pad and Floor mat monitor

**A-BME-03** 



The A-BME-03 Cura1 hardwired monitor is an optional accessory for hardwired floor sensor mat use but is essential for hardwired bed and chair exit sensor pads.

- Pause mode to attend to a patient without triggering alarm.
  - Tamper proof reset option using magnetic caregiver key.
    - Choice of tone, and volume control for alarm.
      - 0, 2s, and 15min. alarm delay settings.
    - Port to connect with nurse call system (NC/NO).
      - Operates on batteries or via AC adaptor.
        - Low battery notification.
    - Input lost notification (pad lost, nurse call, AC lost).
      - Change pad (SafeTime) indication.
- Comes with mounting bracket, hanger strap, soft cover, reset key and batteries.

SAFE-LIFE 5/ 270 LOWER DANDENONG RD MORDIALLOC, VIC 3195

> PH: 03 9588 0833 FAX: 03 9588 0933