

RFID SECURITY

INCREASE SAFETY FOR YOUR PATIENTS AND RESIDENTS

Patient Wander Alert System

Protrac iD Patient Wandering Alert Systems detect wandering patients wearing Protrac iD RFID tags when they approach monitored doors or areas both inside and outside.

Depending on the settings and the mode of operation selected the PWAS triggers alarms, locks or open doors, alerts staff and communicates with other systems such as Nurse Call, Access Control, Radio Paging, DECT Phones, Protrac iD AssetTrac or BabyGuard[™] Software.

Protrac iD Patient Wandering Alert System is available as a 2R5002 or 2R5502 and the 2R2800PWAS is available where multiple doors (up to 4) require monitoring in the same facility.

Protrac iD tags are individually identified and transmit a "heart beat" radio signal continuously with a battery life exceeding 3 years.

Benefits of the Protrac iD Patient Wander Alert System:

- Give more freedom to patients/residents
- Integrated into most new & existing Nurse Call Systems
- Indoor & Outdoor alert monitoring
- Up to 150m detection range
- No ugly portals
- Wireless
- Door lock
- Staff override



2R5002PW



2R5502PW



1RN350 Patch Antenna



2T7000M Bracelet Tag



Distributed by:

🐼 (E 🕑 FC

Designed and Manufactured in Australia. All rights reserved. Specifications subject to change without notice. www.protracid.com.au