

CORDLESS FLOOR MAT

A-MAXI-MAT-CL



The Cordless Floor mat is to be used with the A-CL433-TX (Cordless Pad Transmitter), in order to detect unassisted movement from a bed.

The Cordless floor mat is made from a heavy-duty, non-slip material top (Polyurethane) and the bottom made from PVC.

The cordless mat has a discreet pocket to insert the transmitter (A-CL433-TX) into to avoid damage.

The colour is a charcoal grey/ Black colour with a "Yoga Top" surface.

When a resident removes themselves from bed and steps onto the floor mat the transmitter will send an alarm to the monitor (A-CBM-02).

The Cordless floor mat is available in kits to suit home care and residential age care applications.

Technical Specifications:

- Dimensions 1220 x 610mm
- 1 Year life expectancy
- Recommended weight range for reliable operation: 30kg to 200kg

Smart Caller Safe Life
5/ 270 Lower Dandenong Rd Mordialloc Vic
3195
PH: +61 3 9588 0833
E: sales@safelife.com.au